



## Aerodyne Mini-TILDAS CO<sub>2</sub> Isotope Monitor

*Unprecedented precision and time response for  $\delta^{13}\text{C}$  and  $\delta^{18}\text{O}$  isotopic ratios in a compact, transportable package.*




### Features:

- $< 0.10$  ‰ precision for  $\delta^{13}\text{C}$ ,  $\delta^{18}\text{O}$  in 1 s
- $< 0.03$  ‰ precision for  $\delta^{13}\text{C}$ ,  $\delta^{18}\text{O}$  in 100 s
- Fast time response (10 Hz)
- Direct measurement of CO<sub>2</sub> isotopes in air without sample processing
- Dual laser package allows simultaneous measurement of water isotopes or 12C17O16O or “clumped” CO<sub>2</sub>: 13C18O16O.
- Repeatability exceeding 0.01 per mil for  $\delta^{13}\text{C}$  and  $\delta^{18}\text{O}$  for a
- 30 minute measurement including balanced working reference measurements.

### TILDAS TECHNOLOGY

Aerodyne instruments use tunable infrared laser direct absorption spectroscopy (TILDAS) at mid-IR wavelengths to probe molecules at their strongest “finger-print” transition frequencies. We further enhance sensitivity by employing a patented multi-pass broad-band absorption cell that provides optical path lengths up to 76 m. Direct absorption spectroscopy allows for fast ( $< 1$  sec) absolute trace gas concentrations without need for elaborate calibration procedures. Moreover, TILDAS instruments are relatively free of measurement interference from other molecular species, enabling extremely specific detection.



### Rugged, field-ready instruments

*Direct absorption spectroscopy allows for highly specific and accurate gas detection*

*Mid-IR detection enables maximum measurement sensitivity*

### APPLICATIONS

- Determination of atmospheric sources, sinks, and transport through CO<sub>2</sub> isotopic ratios.
- Biosphere exchange.
- Analysis of CO<sub>2</sub> samples derived from marine carbonate.
- Laboratory measurements of discrete samples.
- Carbon capture and sequestration monitoring.
- Breath analysis.

### AERODYNE CO<sub>2</sub> ISOTOPE ADVANTAGES

- Measurement precision comparable to much larger and more expensive IRMS instruments.
- Time response up to 10 Hz enables eddy covariance studies.
- Powerful TDLWintel software provides flexible instrument control, and real-time data analysis.
- Valve control capable of complex scheduling, and automatic backgrounds and calibrations.
- 19” rack mountable for easy installation.
- Turn-key design allows unattended operation in remote field sites.

# Mini-TILDAS CO<sub>2</sub> Isotope Monitor

## SPECIFICATIONS

Isotope ratio precision (1 $\sigma$ )\*

	1 sec	100 sec
$\delta^{13}\text{C}$	0.1 ‰	0.03 ‰
$\delta^{18}\text{O}$	0.1 ‰	0.03 ‰

\*Accuracy before calibration is typically 20‰

Species Precision (1 $\sigma$  @ 400 ppm CO<sub>2</sub>)

	1 sec	100 sec
CO <sub>2</sub>	25 ppb	10 ppb

Time Response

1-10 Hz data rate  
0.05 s minimum Rise/Fall time (1/e)  
(depends on vacuum pump)

Dynamic Range (air)

	min	max
CO <sub>2</sub>	0 ppm	5,000 ppm


Related Instruments

Dual laser isotope monitor for  $\delta^{13}\text{C}$ ,  $\delta^{18}\text{O}$  and  $\Delta^{17}\text{O}$  of CO<sub>2</sub>

Single laser isotope monitor for  $\delta^{18}\text{O}$  and  $\Delta^{17}\text{O}$  of CO<sub>2</sub>

Dual laser monitor for CO<sub>2</sub> ( $\delta^{13}\text{C}$ ,  $\delta^{18}\text{O}$ ) and water ( $\delta^{18}\text{O}$ ,  $\delta\text{D}$ ) isotopes

Experimental spectrum acquired at 1 Hz



Installation

19" rack mountable or benchtop  
Flushing the optics with CO<sub>2</sub>-free gas is recommended

Instrument Operations

Operating temperature: 10 to 35 °C  
Sample flow rate: 0 to 20 slpm

Instrument Components

Core instrument  
Thermoelectric chiller  
Keyboard, mouse, and monitor  
Vacuum pump (customer specified)  
Inlet sampling system (customizable)

Data Outputs

RS-232, USB drive, ethernet

Size, Weight, Power

Dimensions: 440 mm x 660 mm x 6U (267mm) (W x D x H)  
Weight: 35 kg (core instrument) + 15 kg (chiller) + pump weight  
Electrical Power: 250 W, 120/240 V, 50/60 Hz (without pump)

*Aerodyne specializes in collaboration and custom design. Please contact us if you would like to discuss additional measurement options and applications.*

## REFERENCES

Wehr, R., Munger, J.W., McManus, J.B., Nelson, D.D., Zahniser, M.S., Davidson, E.A., Wofsy, S.C. and Saleska, S.R., Seasonality of temperate forest photosynthesis and daytime respiration, *Nature*, 534(7609), 680-683, 2016.

Sakai, S., Matsuda, S., Hikida, T., Shimono, A., McManus, J.B., Zahniser, M., Nelson, D., Dettman, D.L., Yang, D. and Ohkouchi, N., High-Precision Simultaneous 18O/16O, 13C/12C, and 17O/16O Analyses for Microgram Quantities of CaCO<sub>3</sub> by Tunable Infrared Laser Absorption Spectroscopy. *Analytical chemistry*, 89(21), 11846-11852, 2017.

McManus, J. Barry, David D. Nelson, and Mark S. Zahniser, Design and performance of a dual-laser instrument for multiple isotopologues of carbon dioxide and water, *Optics express* 23, 5, 6569-6586, 2015.

Wehr, R., et al., Long-term eddy covariance measurements of the isotopic composition of the ecosystem-atmosphere exchange of CO<sub>2</sub> in a temperate forest, *Agricultural and forest meteorology*, 181, 69-84, 2013.