



CO₂ Isotope Monitor for $\delta^{13}\text{C}$, $\delta^{18}\text{O}$ and $\Delta^{17}\text{O}$

Direct Spectroscopic Measurement of $\Delta^{17}\text{O} - \text{CO}_2$ with No Chemical Processing.




Features:

- < 0.10 ‰ precision for $\delta^{13}\text{C}$ and $\delta^{18}\text{O}$ in 1s.
- < 0.15 ‰ precision for $\delta^{17}\text{O}$ in 1s.
- Direct measurement of CO₂ isotopes in air without sample processing.
- Repeatability exceeding 0.015 ‰ for $\delta^{13}\text{C}$, $\delta^{17}\text{O}$ and $\delta^{18}\text{O}$ for a 30-minute measurement including balanced working reference measurements.
- Suitable for CO₂ samples derived from carbonate via acid digestion.

TILDAS TECHNOLOGY

Aerodyne instruments use tunable infrared laser direct absorption spectroscopy (TILDAS) at mid-IR wavelengths to probe molecules at their strongest “finger-print” transition frequencies. We further enhance sensitivity by employing a patented multi-pass broad-band absorption cell that provides optical path lengths up to 400 m. Direct absorption spectroscopy allows for fast (<1 sec) absolute trace gas concentrations without need for elaborate calibration procedures. Moreover, TILDAS instruments are relatively free of measurement interference from other molecular species, enabling extremely specific detection.



Rugged, field-ready instruments

Direct absorption spectroscopy allows for highly specific and accurate gas detection

Mid-IR detection enables maximum measurement sensitivity

APPLICATIONS

- Determination of atmospheric sources, sinks, and transport through CO₂ isotopic ratios.
- Biosphere exchange.
- Analysis of CO₂ samples derived from marine carbonate.
- Laboratory measurements of discrete samples.
- Carbon capture and sequestration monitoring.
- Breath analysis.

AERODYNE CO₂ ISOTOPE ADVANTAGES

- Measurement precision comparable to much larger and more expensive IRMS instruments.
- Time response up to 10 Hz enables eddy covariance studies.
- Powerful TDLWintel software provides flexible instrument control, and real-time data analysis.
- Valve control capable of complex scheduling and automatic background and calibrations.
- Optional automated sample handling systems.
- Turn-key design allows unattended operation.

CO₂ Isotope Monitor for $\delta^{13}\text{C}$, $\delta^{18}\text{O}$ and $\Delta^{17}\text{O}$

SPECIFICATIONS

Continuous Measurement Specifications – High Speed

	CO ₂	$\delta^{13}\text{CO}_2$	$\delta\text{CO}^{18}\text{O}$	$\delta\text{CO}^{17}\text{O}$
0.1 second	0.15 ppm	0.3‰	0.3‰	0.5‰
1 second	0.05 ppm	0.1‰	0.1‰	0.15‰
60 second	0.015ppm	0.03‰	0.03‰	0.05‰

Note: These measurements are not referenced to a working reference. This configuration supports 10 Hz eddy covariance measurements with a modest sized pump (120 lpm) and a flow rate of 6 slpm.

Continuous Measurement Specifications – High Precision

	CO ₂	$\delta^{13}\text{CO}_2$	$\delta\text{CO}^{18}\text{O}$	$\delta\text{CO}^{17}\text{O}$
2 min measurement	0.02 ppm	0.03‰	0.03‰	0.05‰
20 min measurement	0.01 ppm	0.01‰	0.01‰	0.015‰

Note: These measurements are normalized to a working reference and the time to do so is included in the quoted measurement time. The working reference has a mixing ratio similar to the sample. The flow rate is 0.6 slpm.

Discrete Sample Specifications – High Precision


	CO ₂	$\delta^{13}\text{CO}_2$	$\delta\text{CO}^{18}\text{O}$	$\delta\text{CO}^{17}\text{O}$
One air sample (~15 ml of air or 0.25 $\mu\text{moles CO}_2$) 3 min measurement	0.02 ppm	0.03‰	0.03‰	0.05‰
10 Air samples (~150 ml of air or 2.5 $\mu\text{moles CO}_2$) 30 min measurement	0.01 ppm	0.01‰	0.01‰	0.015‰

Note: These measurements are normalized to a working reference and the time to do so is included in the quoted measurement time. The working reference has a mixing ratio similar to the sample.

Related Instruments

Single laser isotope monitor for $\delta^{13}\text{C}$ and $\delta^{18}\text{O}$ of CO₂
Single laser isotope monitor for $\delta^{18}\text{O}$ and $\Delta^{17}\text{O}$ of CO₂
Dual laser monitor for CO₂ ($\delta^{13}\text{C}$, $\delta^{18}\text{O}$) and water ($\delta^{18}\text{O}$, δD) isotopes

Experimental Spectrum Acquired at 1 Hz



Installation

19" rack mountable or benchtop
Flushing the optics with CO₂-free gas is recommended

Instrument Operations

Operating temperature: 10 to 35 °C
Sample flow rate: 0 to 20 slpm

Instrument Components

Core instrument
Thermoelectric chiller
Keyboard, mouse, and monitor
Vacuum pump (customer specified)
Inlet sampling system (customizable)

Data Outputs

RS-232, USB drive, ethernet

Size, Weight, Power

Dimensions: 560 mm x 770 mm x 640 mm (W x D x H)
Weight: 75 kg
Electrical Power: 250-500 W, 120/240 V, 50/60 Hz (without pump)

Aerodyne specializes in collaboration and custom design. Please contact us if you would like to discuss additional measurement options and applications.

REFERENCES

Wehr, R., Munger, J.W., McManus, J.B., Nelson, D.D., Zahniser, M.S., Davidson, E.A., Wofsy, S.C. and Saleska, S.R., Seasonality of temperate forest photosynthesis and daytime respiration, *Nature*, 534(7609), 680-683, 2016.

Sakai, S., Matsuda, S., Hikida, T., Shimono, A., McManus, J.B., Zahniser, M., Nelson, D., Dettman, D.L., Yang, D. and Ohkouchi, N., High- Precision Simultaneous 18O/16O, 13C/12C, and 17O/16O Analyses for Microgram Quantities of CaCO₃ by Tunable Infrared Laser Absorption Spectroscopy. *Analytical chemistry*, 89(21), 11846-11852, 2017.

McManus, J. Barry, David D. Nelson, and Mark S. Zahniser, Design and performance of a dual-laser instrument for multiple isotopologues of carbon dioxide and water, *Optics express* 23, 5, 6569-6586, 2015.

Wehr, R., et al., Long-term eddy covariance measurements of the isotopic composition of the ecosystem-atmosphere exchange of CO₂ in a temperate forest, *Agricultural and forest meteorology* 181, 69-84, 2013.